In The Name of God







Monitoring ECG Signal Instrumentation Project

Group Members:

Maral Mordad


Sina Rashidizade

Mahya Haghgoo

Mohadese Ghsemmehrabi


## Steps






# Introduction What is ECG?


An electrocardiogram records the electrical signals in your heart.



# ECG How can be measured?



# ECG The Main Pick






## ECG Disease Detection

- Arrhythmias where the heart beats too slowly, too quickly, or irregularly
- Coronary heart disease where the heart's blood supply is blocked or interrupted by a build-up of fatty substances
- Heart attacks where the supply of blood to the heart is suddenly blocked

# ECG Features: Heartbeat / HRV



### ECG Noises

- $\mathbb{W}-$  Noise in ECG data

**High-Frequency Noise**: Power line Noise, White Gaussian Noise, Eletcromyogram/Motion Noise

Low-Frequency Noise: Baseline drift, Electrode contact loss

IIR Notch Filters FIR Filters



#### Data Extraction

\_\_\_\_\_ Data Base : www.archive.physionet.org


\_/<sub>M</sub> Data Information :

**ECG Signal** 

MIT-BIH Normal Sinus Rhythm Database(nsrdb)


No. of samples: 600


Sampling Frequency: 128




### Processing

Pan-tompkins algorithm






#### Circuites



#### Proteous Environment




## Raspberry

```
import time
      import serial
      import requests
      ser = serial.Serial(port='/dev/ttyS2', baudrate=9600, timeout=1)
    ⊟while 1:
          x=ser.readline()
          x = x.decode('ascii')
          if x == '':
10
11
              continue
12
          x = x.strip('\r\n')
13
14
          try:
              eq = x.find('=')
15
              num = int(x[eq+1:])
16
17
          except:
18
              continue
19
          #print(num)
20
21
          requests.get('http://127.0.0.1/send?heart={}&ecg={}'.format(3, num))
22
23
```

## Webserver Before Starting





### Webserver



3 192.168.43.12/ui








```
sampled data | Arduino 1.8.13
                                                                                                                                                     File Edit Sketch Tools Help
                                                                                                                                                         Ø-
 sampled data
 4 //float data[1280] = {-37, -39, -39, -41, -27, -21, -27, -21, -27, -27, -31, -27, -25, -27, -37, -33, -25, -29, -23, -25, -31, -29, -33, -29, -27, -31, -31, -3
 5 int i = 0:
 6 const int datasize = 600;
 7 int data[datasize] = {-33, -31, -39, -41, -37, -31, -27, -19, -15, -13, -13, -25, -25, -25, -25, -25, -23, -25, -23, -27, -27, -27, -25, -35, -29, -25, -29, -27, -3
9 void setup() {
10 Serial.begin (9600);
11
12 }
13
14 void loop() {
15
    char buffer [15];
16
17
18
    sprintf (buffer. "data=%d". data[i++]);
19
20 Serial.println(buffer);
   if(i>= datasize)
21
22
      i = 0;
23
24
   delay(100);
25
26 }
"B:\\Arduino\\hardware\\tools\\avr/bin/avr-gcc-ar" rcs "C:\\Users\\maral\\AppData\\Local\\Temp\\arduino build 901568\\core\\core.a" "C:\\Users\\maral\\AppData\\Local\\Temp\\
Archiving built core (caching) in: C:\Users\maral\AppData\Local\Temp\arduino cache 549023\core\core arduino avr uno 71f1c74518d18c55e877003a61b70c8c.a
Linking everything together ...
"B:\\Arduino\\hardware\\tools\\avr/bin/avr-qcc" -w -Os -g -flto -fuse-linker-plugin -Wl,--gc-sections -mmcu=atmega328p -o "C:\\Users\\maral\\AppData\\Local\\Temp\\arduino b
"B:\\Arduino\\hardware\\tools\\avr/bin/avr-objcopy" -0 ihex -j .eeprom --set-section-flags=.eeprom=alloc,load --no-change-warnings --change-section-lma .eeprom=0 "C:\\Users
"B:\\Arduino\\hardware\\tools\\avr/bin/avr-objcopy" -O ihex -R .eeprom "C:\\Users\\maral\\AppData\\Local\\Temp\\arduino build 901568/sampled data.ino.elf" "C:\\Users\\maral
```


Arduino Uno on COM5


"B:\\Arduino\\hardware\\tools\\avr/bin/avr-size" -A "C:\\Users\\maral\\AppData\\Local\\Temp\\arduino build 901568/sampled data.ino.elf"

### Mobile





## Thingsboard





## Thingsboard





### Thanks

